**EDU-CHEM INNOVATIONS**

101 Malaga Street  
Royal Palm Beach, FL 33411

Toll Free: (866) 859-9904

Phone: (561) 337-4665

Fax: (561) 828-0563

Email: [sales@edu-chem.com](mailto:sales@edu-chem.com)

Visit our Website @ [www.edu-chem.com](http://www.edu-chem.com)

**Quotation**


To:

Company: JUNIATA COLLEGE  
Attn: TARA FITZSIMMONS  
Email: [fitsit@juniata.edu](mailto:fitsit@juniata.edu)  
Telephone: 814-641-3543

QUOTATION NUMBER: 709140

DATE: APRIL 18, 2007

QUANTITY	ITEM #	DESCRIPTION	UNIT PRICE
1 - 3 UNITS	EDU-SPEC-6405	SINGLE BEAM UV/VISIBLE RANGE SCANNING SPECTROPHOTOMETER (190 - 1100nm) WITHOUT PRINTER	4238.00
4 - 6 UNITS			4068.50
7 - 9 UNITS			3983.75
10 + UNITS			3898.50
1 - 3 UNITS	EDU-SPEC-6412	PRINTER OPTION (1) FACTORY FITTED INTERNAL "CASH REGISTER" TYPE PRINTER	1034.00
4 - 6 UNITS			992.65
7 - 9 UNITS			971.95
10 + UNITS			951.30
1 - 3 UNITS	EDU-SPEC-2828	PRINTER OPTION (2) EXTERNAL SERIAL PRINTER AND CALBE	395.00
4 - 6 UNITS			379.20
7 - 9 UNITS			371.30
10 + UNITS			363.40
1 - 3 UNITS	EDU-SPEC-6442	OPTIONAL EIGHT POSITION MOTORIZED MULTI CELL CHANGER	965.00
4 - 6 UNITS			926.40
7 - 9 UNITS			907.10
10 + UNITS			887.80

**EDU-CHEM INNOVATIONS**

101 Malaga Street  
Royal Palm Beach, FL 33411

Toll Free: (866) 859-9904

Phone: (561) 337-4665

Fax: (561) 828-0563

Email: [sales@edu-chem.com](mailto:sales@edu-chem.com)

Visit our Website @ [www.edu-chem.com](http://www.edu-chem.com)

## Quotation

<b>To:</b>
Company: JUNIATA COLLEGE
Attn: TARA FITZSIMMONS
Email: <a href="mailto:fitsit@juniata.edu">fitsit@juniata.edu</a>
Telephone: 814-641-3543

QUOTATION NUMBER: 701250

DATE: MAY 21, 2007

QUANTITY	ITEM #	DESCRIPTION	UNIT PRICE
1 - 3 UNITS	EDU-SPEC-6310	MODEL 6310 VISIBLE RANGE SCANNING SPECTROPHOTOMETER (320-1000NM) SUPPLIED WITH MAINS LEAD, PACK 100 DISPOSABLE CUVETTES, 10X10MM CELL HOLDER AND INSTRUCTION MANUAL.	3392.00
4 - 6 UNITS			3256.30
7 - 9 UNITS			3188.50
10 + UNITS			3120.65
1 - 3 UNITS	EDU-SPEC-6310PRINT	SPECTROPHOTOMETER PRINTER W/ CABLE	395.00
4 - 6 UNITS			379.20
7 - 9 UNITS			371.30
10 + UNITS			363.40



[Web](#) [Images](#) [Video](#) [News](#) [Maps](#) [Desktop](#) [more »](#)

juniata

Search

[Desktop Preferences](#)  
[Advanced Search](#)**Cached messages**

Message 13 of 13 in conversation

[« Older](#) | [Newer »](#) | [View Entire Thread \(13\)](#) | [Reply](#) | [Reply to all](#) | [Forward](#) | [View in Outlook](#) **Fwd: RE: Fw: tech question about peak simple**

From: &lt;chemchek@aol.com&gt;

To: John (& Annette) Mellor <[john@bucksci.com](mailto:john@bucksci.com)>, Derek (& Molly) Hodgeman <[derek@bucksci.com](mailto:derek@bucksci.com)>, [stephen@bucksci.com](mailto:stephen@bucksci.com), Eric Anderson <[eric@bucksci.com](mailto:eric@bucksci.com)>

Date: Apr 04 2007 - 8:27am

See! I've never had good success with those USB Adapters...

FYI,

DeMento in Agony Mountains

\*\*\*\*\*

-----Original Message-----

From: "Bernauer, Bob"

Subj: RE: Fw: tech question about peak simple

Date: Wed Apr 4, 2007 9:05 am

Size: 12K

To:

Thanks for the response.

It is an issue with an emulator which for some reason I cannot get to work. Easiest is to change computers to one with a com port so that is next.

Thanks again

Bob Bernauer

From: [ChemChek@aol.com](mailto:ChemChek@aol.com) [<mailto:ChemChek@aol.com>]

Sent: Wednesday, April 04, 2007 1:12 AM

To: Bernauer, Bob

Cc: Eric Anderson / Sales &amp; Mktng; Stephen Zakrewski / Inside Sales; John Mellor /

GC-LC Department

Subject: Re: Fw: tech question about peak simple

Bob,

If you have the Manual... this is covered in the beginning of thPeak-Simple Section... and means you do not have the correct Comm-Port selected and the GC / PC cannot communicate.

You have to press several times and then "OK" when the "ACQUISITION" Screen pops up with the count-down Timer to EXIT this mode... then go into "EDIT", "OVERALL" and select the proper Comm-Port. Always SAVE the updated "Control File" whenever you make a change... as these are the parameters that are loaded with the Software.

I am enclosing one of my Mini-Guides for the GC / Peak-simple Operation if you do not have a chance to go through the complete Manual (and please reply to [JOHN@BUCKSCI.COM](mailto:JOHN@BUCKSCI.COM) if you do not have the Manual and he'll get one out to you).

Are you using an ACTUAL Comm-Port or a "Virtual Adapter" on a USB Port... as the latter often does not provide TRUE EMULATION of an RS-232 Serial Port. Maybe your WXP is having a problem if using an Emulator, as well. Call John at BUCK (800-562-5566) if this is the case.

Call me, too, if you still have these set-up problems... as usually we have no difficulties when new Users start up with the proper conditions... and I'm running around the Mid-West to the North-East covering the Natl ACS / Chicago, NSTA / St. Louis, Freshman Gen Chem / U. of Cincinnati, Science-in-Motion Workshop / OSU, Columbus and tomorrow at Juniata College in the middle of the PA Pcono mountains... so will get a Phone Call before my E-Mail for a few days.

Regards,

Dr. Jerry DeMenna, Ph.D.  
General Manager / App-Lab Specialist  
CELL: 800-531-CHEM  
\*\*\*\*\*

----- Original Message -----

From: Bernauer, Bob  
To: [sales@bucksci.com](mailto:sales@bucksci.com)  
Sent: Friday, March 30, 2007 11:37 AM  
Subject: tech question about peak simple

I'm on the learning curve with the peak simple software and the model 910 GC we recently purchased.

Every 15 seconds or so a window pops onto the screen with the heading;

Data Acquisition System  
And it keeps  
Waking up  
Signing on  
Setting modules  
Calibrating

While this occurs I have no input available which is very inconvenient if I have just injected and want to run the program.

Is this normal for the software?

Thank you

Robert Bernauer  
Bellarmine university  
Louisville, ky  
[bbernauer@bellarmine.edu](mailto:bbernauer@bellarmine.edu)

\*\*\*\*\*

See what's free at <http://www.aol.com>.



This e-mail was scanned for viruses using BitDefender

[« Older](#) | [Newer »](#)   [View Entire Thread \(13\)](#)   [Reply](#) | [Reply to all](#) | [Forward](#) | [View in Outlook](#)

[Google Desktop](#) - [Home](#) - [Browse Timeline](#) - [Index Status](#) - [Privacy](#) - [About](#) - ©2007 Google